

STOP!

- Do not attempt to install this safety device without the knowledge of the magnetic field you are working around.
- Any and all tools, screws, etc. are subject to and should be considered to be potentially hazardous not only to the equipment but to the installer, due to the magnetic field. Tools used should be approved for the magnetic field.
- The Safety Stop should always be installed by a qualified installer who has knowledge of the magnetic field and has been properly trained in magnetic safety.
- The Safety Stop should only be installed as part of the door frame.
- The Safety Stop should be installed with 2-1/2" wood screws if the door jamb is of standard wood construction. If door jamb is made of other materials, or if you are unsure of the jamb construction consult with the shielding manufacturer for suggested types of fasteners.
- Do not attempt to install without first reading the entire installation manual.

Installation Instructions

Step 1: The door to the MRI machine should be closed during the installation of the Safety Stop.

Step 2: We recommend installing the Safety Stop in the middle or bottom of the door frame. The instructions below assume you are installing in the center. See drawing on reverse side of this page.


Step 3: To install top bracket, first find the center of the door jamb top to bottom. Next, measure up from the center of the door jamb 6-1/4". This is the bottom of the top bracket. Install top bracket with 2-1/2" wood screws. Make sure the sign is 100% vertically and horizontally level!

Step 4: Place sign onto top bracket.

Step 5: Place bottom bracket on top of the top pins. Make sure to keep spacing between top and bottom bracket 12-1/2" to 12-3/4". Install top hinge bracket to the door jamb, allowing for the swing of the sign. It is recommended to test the swing of the sign before installing all screws.

Step 6: Before opening MRI door, check to make sure that all tools, and any magnetic objects are removed from the area.

SCALE = 1/2" = 1'



SHEET 1 OF 1	PROJECT NO.	SAFETY STOP SIGN		DATE DRAWN: MAY 25 2007
		MOUNTING SPECIFICATIONS		DRAWN BY: MEG
THESE PRINTS ARE INTENDED FOR THE EXCLUSIVE USE BY PRECISION MILLWORK & PLASTICS, INC.				PRODUCTION MANAGER APPROVAL: OK
				PROJECT MANAGER APPROVAL: OK