



Examine Real Quality & Value™

Laminate

Cleaning Tips

To clean the surface, use a damp cloth or sponge and a mild soap or detergent. Difficult stains such as coffee or tea can be removed using a mild household cleaner and baking soda; mixing to achieve a paste consistency. Use a stiff nylon bristle brush, scrubbing (approx. 15-20 strokes) the affected area. Do not scrub so as to mar (damage, scratch) the surface finish.

Stubborn stains that resist any of the above cleaning methods may require the use of undiluted household bleach or nail polish remover. Use a cotton ball saturated with bleach or nail polish remover (acetone), gently rub the stain for up to two minutes. Rinse thoroughly with warm water and wipe dry using a soft cloth. This step may be repeated if the stain appears to be going away and the color of the laminate has not been affected.

WARNING: Prolonged exposure of the laminate surface to bleach will cause discoloration. Always rinse laminate surfaces after cleaning! Failure to rinse after cleaning can cause damage; even if a small amount of cleaning solution remains on the surface. A dry residue may be invisible; however, moisture from cups or drinks can reactivate it, and result in permanently etched scars or stains over time.

Recommended Household Cleaners:

Formula 409®

Lestoil®

Dawn®

Lysol® Brand Disinfectant Basin/Tub/Tile Cleaner

Glass Plus®

Mr. Clean®

Dow Bathroom Cleaner with Scrubbing Bubbles™

TOP JOB®

Fantastik®

Clorox®*

Windex®

* Prolonged exposure can mar the surface

A Few Notes of Caution

Acidic or **abrasive** cleaners can damage laminate surfaces; do not use them. Drain cleaners containing lye will permanently damage the laminate surface. If you spill a drain cleaner, wipe it up immediately and rinse several times with water.

Hair, textile and food **dyes** can cause permanent stains. If dye should happen to spill, wipe it up immediately with dishwashing detergent or an all-purpose cleaner.

Wipe spills away promptly and rinse several times with water.

Rust removers contain harsh chemicals, which will quickly cause permanent damage. If a spill occurs, wipe off all residues immediately; wash thoroughly with soapy water and rinse several times.

Steel wool and other abrasive pads will damage the laminate face. Do not use them for cleaning and don't store steel wool pads on your countertop; the metal can rust and leave stains.

Toilet bowl cleaners contain harsh chemicals that can cause permanent damage. If spills occur, wipe up immediately, wash surface with soapy water and rinse several times.

The cleaners listed below can cause damage to the surface of laminate:

Chemical Ingredient	Synonymous Names
Hydrochloric Acid	Muriatic Acid Hydrogen Chloride
Sulfuric Acid	Oleic Acid Oil of Vitriol Oleum
Hydrofluoric Acid	Rust Remover
Phosphoric Acid	Rust Remover
Sodium Hydroxide	Caustic Soda Caustic Lye Soda Lye
Pumice (abrasive)	

Remember, sharp objects can damage the surface of laminate; marring its beauty and lowering wear and stain resistance. Although high-pressure laminates are somewhat resistant to scratch and marring, the surface can be damaged, even under normal use and is not covered under warranty.

NOTE: Clinton Industries is not responsible for damage to the laminate surface as a result of proper or improper attempts at cleaning the surface of the laminate. Any damage to the surface or coloration of the laminate is the sole responsibility of the purchaser of the product. Clinton Industries, Inc. is not liable to replace or repair damaged laminate as a result of attempts to clean the surface.