

# Instruction Sheet

## Computer Security Cart #350717 - #350718


- 1) The basket should contain:  
(4) 10-32 x 3/8 Countersunk Head Screws - (4 sets) Dual Lock Mounting Tape (1) Poly Tray - (4) 4 x 12 screws
- 2) Cut the straps holding the basket to the center column of the cart and remove the basket.
- 3) Remove and unwrap the Monitor Plate Head Assembly.
- 4) With a second person, grasp the cart by the base and ring handle, lift it from the carton, and set it upright on the floor.
- 5) Place the basket hook into the receiver loop on the column and the poly tray into one of the quadrants of the ring handle.
- 6) Use the 4 provided flat head screws to attach the Monitor Plate Head Assembly to the battery housing.
- 7) Wipe the surface of the Keyboard Tray with isopropyl alcohol.
- 8) Remove the backing from one side of each of the Dual Lock tape sets. Attach each to an appropriate location on the bottom of the keyboard.
- 9) Remove the backing from the other side of the Dual Lock and position the keyboard on the Keyboard Tray. Press it firmly into place.
- 10) Remove the four screws that hold the access plate in place.
- 11) Insert all cables and wires and organize in the cable management housing.
- 12) Re-install the access plate. Firmly tighten screws in place.
- 13) Use M4 screws from monitor to secure the monitor to the monitor plate.

### Height Adjustment

To raise stand, press down on foot ring. Note: When not in use for a period of time, use both hands with equal pressure and press down on ring handle slightly while stepping on foot ring. To lower stand, use both hands with equal pressure and press down on ring handle while stepping on foot ring.

### Keyboard Adjustment


Adjust keyboard tray as needed by loosening the knobs provided




MANUFACTURERS OF **IBeam**® PRODUCTS  
**OMNIMED** INC.®

# Instruction Sheet

## CPU Holder



Use the provided security screw wrench to loosen the Top Cover off the CPU Holder



Loosen the CPU Slide screw and adjust the width of the slides to match your PC. Loosen the screws holding the side slide in place. Place your CPU on the base and position the slides closely around it. Remove the CPU and tighten the slide screws. Replace the CPU and secure it in by retightening the Top Cover.

Your Transport Stand is now fully assembled and ready for use.

MANUFACTURERS OF **iBEAM**® PRODUCTS  
**OMNIMED** INC.®