...going one step further

SONOtrain™

1019634, 1019635, 1019636, 1019637, 1019638

P120 SONOtrain™ Vein model

Ultrasound block with three blood vessels with diameters of 4, 8 and 15 mm with adjustable fluid flow. Punctures and injections possible, and self-closing injection channels.

Replacement parts:

1019652 P120: Replacement block for vein model

Cleaning instructions:

Ultrasound scan gel residue can be removed with a damp cloth or under running water. After drying the model, sprinkle on the talcum powder provided. Rub it over the entire surface and then put the block away in the storage box.

Technical data:

Dimensions: Block (225 x 180 x 50 mm), storage box (360 x 360 x 125 mm)

Weight: Block (approx. 1600 g), storage box with contents (approx. 3200 g)

Operating conditions:

Operating temperature 0°C to +40°C (32°F - 104°F)

Storage temperature -10°C to +60°C (14°F - 140°F) (only store in the box provided!)

The simulator contains no substances indicated in the REACH regulation EC no. 1907/2006

Note:

The accessories for building the model are in the space under the block. This space can also be used as a container to catch any fluids used.

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P121 SONOtrain™ Foreign body model

English

Ultrasound block with 6 different foreign bodies for insertion: Pellets, shard of glass, nail, projectile, loose chippings and splinters of wood can be clearly detected and recognized.

Replacement parts:

1019651 P121: Replacement block for foreign body model

Cleaning instructions:

Ultrasound scan gel residue can be removed with a damp cloth or under running water.

After drying the model, sprinkle on the talcum powder provided. Rub it over the entire surface and then put the block away in the storage box.

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Weight: Block (approx. 1800 g), storage box with contents (approx. 3000 g)

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Breast model with cysts

English

Realistic reproduction of a breast made from ultrasound material with two cysts for realistic punctures and aspirations.

Replacement parts:

1019649 P124: Replacement breast with cysts

Cleaning instructions:

Ultrasound scan gel residue can be removed with a damp cloth or under running water. After drying the model, sprinkle on the talcum powder provided. Rub it over the entire surface and then put the breast away in the storage box.

Technical data:

Dimensions: Breast (170 x 130 x 55 mm), storage box (360 x 360 x 125 mm)

Weight: Breast (approx. 500 g), storage box with contents (approx. 1700 g)

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P125 SONOtrain™ Breast model with tumours

Realistic reproduction of a breast made from ultrasound material with three tumours for realistic punctures and biopsies.

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After drying the model, sprinkle on the talcum powder provided. Rub it over the entire surface and then put the breast away in the storage box.

Technical data:

Dimensions: Breast (190 x 140 x 60 mm), storage box (360 x 360 x 125 mm)

Weight: Breast (approx. 700 g), storage box with contents (approx. 1900 g)

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© Copyright 2015 for instruction manual and design of product:
3B Scientific GmbH, Germany