## 1013116 | P58/1

[www.3bscientific.com](http://www.3bscientific.com)

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# PROSTATE EXAMINATION MODEL

## Introduction

You have purchased a prostate examination model, which can be used to demonstrate and practice the digital rectal examination (DRE). The examination model has 4 different prostates, including their respective sonographic images, which can be examined and reviewed in each case, allowing the training and improvement of standards within the framework of professional screening in specialist medical training and everyday clinical practice.

## Delivery contents

Prostate examination model

Figure 1: Prostate examination model

## Functions

The prostate examination model is used to learn and practice the digital rectal examination (DRE). The index finger is inserted into the rectum, facilitating palpation of the rear part of the prostate, where the majority of prostate carcinomas are found. The shape, position and consistency of the prostate are assessed, allowing irregularities, induration and nodes on the otherwise smooth surface of the prostate to be felt.

The model presents four different prostates including transrectal sonographic imaging (transrectal ultrasound (TRUS)):


Figure 2: Disease entities

### 1. Normal prostate (Figure 2a):

The normal healthy prostate is chestnut-sized, has a sharply defined edge and a smooth surface. Two side lobes with a shallow median furrow can be felt. The consistency is firm-elastic (rubber-like).

### 2. Advanced prostate carcinoma (Figure 2b):

The prostate is hardened. The appearance is characterized by nodes and bumpy irregularities. The prostate cannot be clearly delineated.

### 3. Prostate carcinoma at an early stage (Figure 2c):

A small local induration is palpable. The prostate is clearly delineated and not enlarged.

### 4. Benign prostatic hyperplasia (BPH) (Figure 2d):

The prostate is enlarged and has a firm-elastic consistency. It is clearly delineated and symmetrical.

**NOTE:**

The sonographic images for each prostate are displayed opposite the examined prostate. This means that when the palpation is performed correctly (side position or bent forward, see Figure 3), the associated sonographic image is not in the user's field of vision. This facilitates subsequent result checking and self-monitoring (please note numbering). Please also note also the different perspective of the sonographic imaging compared to palpation (dorsal position for TRUS, see Figure 4).

Figure 3: Palpation with alternative positions

Figure 4: Transrectal ultrasound (TRUS)

**Cleaning and care**

The surfaces of the simulator can be cleaned with a damp cloth. Please ensure that the model is completely dry before storing.

**NOTE:**

Do not use any cleaning products containing solvents, because the surface could be damaged. Please avoid labeling and marking the surface, as this can cause permanent discoloration.

**Technical data**

Dimensions: 9,5 x 9,5 cm  
 Weight: Approx. 0.13 kg (net), approx. 0.28 kg (gross)

**Operating conditions**

Operating temp: -10 °C to +40 °C  
 Storage temp: -20 °C to +60 °C

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© Copyright 2016 for instruction manual and design of product: 3B Scientific GmbH, Germany