

CPR Prompt®+ Powered by Heartisense™ Training & Practice Manikins

Product Name	CPR Prompt®+ Powered by Heartisense™ Training & Practice Manikins
Addresses Training Needs/Curricula	<ul style="list-style-type: none"> • Can be used to train anyone to administer proper CPR technique • Compliant with current American Heart Association guidelines
Configuration	<ul style="list-style-type: none"> • Blue/Tan manikin • Lung bags with insertion tool • Add-on kit • Instruction manual
Product Code/SKU	<p>Adult CPR Prompt®+ (Blue LF06001 A/Tan LF06003 A): Ship weight 7 lbs.</p> <p>5-pk Adult CPR Prompt®+ (Blue LF06100 A/ Tan LF06102 A): Ship weight 29 lbs.</p> <p>Infant CPR Prompt®+ (Blue LF06002 A/Tan LF06012 A): Ship weight 4.75 lbs.</p> <p>5-pk Infant CPR Prompt®+ (Blue LF06050 A/ Tan LF06051 A): Ship weight 28 lbs.</p> <p>Adult CPR Prompt®+ Add-On Kit LF06800: Ship weight 0.25 lbs.</p> <p>5-pk Adult CPR Prompt®+ Add-On Kit LF06801: Ship weight 1.25 lbs.</p> <p>Infant CPR Prompt®+ Add-On Kit (Blue LF06922 A/ Tan LF06931 A): Ship weight 1.25 lbs.</p> <p>5-pk Infant CPR Prompt®+ Add-On Kit (Blue LF06922 B/Tan LF06931 B): Ship weight 7 lbs.</p>



Be
READY

Ready for the unexpected. Ready to save lives.



Powered by heartisense™

Training & Practice Manikins

Professional quality in a cost-effective platform for learning the essentials of basic life support

- American Heart Association compliant
- Get real-time feedback on compressions/ventilations
- Rugged, reliable and durable
- 100% made in the USA
- 5-year warranty on manikin
- 1-year warranty on electronics




CPR Add-On Kit powered by Heartisense™ App
Easily upgrades an existing CPR Prompt® manikin to a CPR Prompt®+.



CPR Prompt®+ Powered by Heartisense™

A professionally made platform that enables effective instruction on the essentials of basic life support including CPR, mouth-to-mouth resuscitation and abdominal thrusts (Heimlich maneuver)

Realistic anatomical markings to learn correct hand placement



Adult model (LF06001 A)

NEW Infant model (LF06002 A)



CPR Add-On Kit powered by Heartisense™ App*

Links up with six manikins at once. Instructors can monitor students using iPhone®, iPad® or Android™ tablets and smartphones with the option to export/print results.

*iPhone®, iPad® and Android™ tablets/smartphones not included. These manikins are intended to be used as an instructional aide in the teaching and demonstrating of live rescue training. Actual products may vary slightly from photos. Nasco Healthcare reserves the right to change product colors, materials, supplies or functions as needed.