






FusionONE Series Exam Table Comparison

All the below exam table models feature ADA low access height and 500-lb weight capacity

Manual Back Exam Tables	Power Hi-Lo	 <p>FusionONE™</p> <ul style="list-style-type: none"> › Power Hi-Lo & Manual Back › Foot control › Front storage drawer › Reversible side storage drawer › Drain Pan <p>Model: 3001</p>	 <p>FusionONE™</p> <ul style="list-style-type: none"> › Power Hi-Lo & Manual Back › Foot control › Front storage drawer › Reversible side storage drawer › Drain Pan <p>Model: 3002</p>	 <p>FusionONE™</p> <ul style="list-style-type: none"> › Power Hi-Lo & Manual Back › Foot control › Front storage drawer › Reversible side storage drawer › Drain pan › Stirrups <p>Model: 3003</p>
		<p>PLUS</p> <ul style="list-style-type: none"> › Stirrups 	<p>PLUS</p> <ul style="list-style-type: none"> › Drawer warmer › Pelvic tilt › Electrical receptacle 	
Power Back Exam Table	Power Hi-Lo	 <p>FusionONE+™</p> <ul style="list-style-type: none"> › Power Hi-Lo & Power Backrest › Foot control › Front storage drawer › Reversible side storage drawer › Drain Pan <p>Model: 3502</p>	 <p>FusionONE+™</p> <ul style="list-style-type: none"> › Power Hi-Lo & Power Backrest › Foot control › Front storage drawer › Reversible side storage drawer › Drain Pan › Stirrups <p>Model: 3503</p>	
		<p>PLUS</p> <ul style="list-style-type: none"> › Stirrups 	<p>PLUS</p> <ul style="list-style-type: none"> › Drawer warmer › Pelvic tilt › Electrical receptacle 	

Comparison chart →



Manual Back Exam Tables

Power Hi-Lo

Power Back Exam Tables

Power Hi-Lo

	3001	3002	3003	3502	3503
Efficiencies					
Power hi-lo with power back				<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Power hi-lo with manual back	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
Foot control	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Pelvic tilt (7°)			<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>
Plastic drain pan	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Large front storage drawer	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Large, reversible side storage drawer	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Foldaway stirrups with four lateral positions, adjustable length		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Infection Prevention					
Removable, seamless upholstery with PreFixx® Protective Finish	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Easy-clean powder coated enamel surfaces	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Patient Experience					
Drawer warmer			<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>
Comfort front leg cushion	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Accessibility					
500lb (226.8kg) patient weight capacity	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Low-access ADA height 17"-19" (43.2cm-48.3cm)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Electrical					
Hospital grade electrical receptacle with child safety cover			<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>
UL60601-1 and CAN/CSA-22.2 No. 601.1 classification	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Quality					
3-year limited warranty	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Durable galvanized steel construction	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>