

# LEAD APRON RACK

INSTRUCTIONS OF USE

## USE

---

Hang leaded apron on hanger (or peg) depending on model.  
respect the limit of one apron per peg/hanger.  
Do NOT overload.

## WARRANTY

---

Phillips Safety Lead Apron Racks have one year warranty against manufacturer defects.  
Attention: Does NOT include normal wear and tear or overloading of hanging points.

## CLEANING

---

Use a mild detergent or cleaner to remove any dirt or debris.  
Do NOT use abrasives.  
Water temperature should be at least 110 degrees.

## DISINFECTING

---

Lead apron racks can be disinfected with most standard hospital grade disinfectants containing:  
QATs, IPA, or Hydrogen Peroxide.