

SureFrame™ PCR Compatibility Chart

UM UNIVERSAL MEDICAL	P9096-HSC	P9096-HSW	P9096-FSC	P9096-FSW	P9384-FSC	P9384-FSW
ABI						
3100, 3130XL, 3500, 3500XL, 3730, 3730XL					•	•
GeneAmp 2700, 2720	•	•				
GeneAmp 9600	•	•				
GeneAmp 9700	•	•			•	•
MiniAmp	•	•				
ProFlex	•	•				
Proflex, Multiblock system					•	•
QuantStudio 5, 6, 7, 12 K Real-time PCR System	•	•			•	•
SimpliAmp	•	•				
Veriti 384-well					•	•
Veriti 96-well	•	•				
ViiA7, 7900 HT FAST					•	•
Agilent						
AriaMx qPCR System			•	•		
Mx3000P, Mx3005P	•	•				
Robocycler Gradient	•	•	•	•		
Surecycler 8800	•	•			•	•
Analytica Jena						
Trobot 96, SpeedCycler2 SP, SPR						
Alpha SC						
Flexcycler2, T1 Thermocycler, T Gradient, T One, T Advanced, T Professional series	•	•	•	•	•	•
qTOWER3 (G, touch), Topical	•	•			•	•
Bioer						
Gene Explorer	•	•				
Gene Touch 384					•	•
Gene Touch 96	•	•				
Biometra						
Biometra T Advanced					•	•
Biometra T Advanced 96, 96 G, 96 S, 96 SG	•	•	•	•		
Biometra T One	•	•	•	•		
Biometra T Professional series	•	•	•	•		
Biometra Uno	•	•	•	•		
BioRad / MJ Research						
C1000, C1000 Touch, S1000	•	•	•	•	•	•
CFX Connect, CFX96, CFX96 Touch			•	•		
CFX384, CFX384 Touch					•	•
CFX96 Touch Deep Well	•	•	•	•		
Chromo4	•	•	•	•		
DNA Engine Opticon and Opticon 2			•	•		
DNA Engine, DNA Engine Tetrad, Tetrad 2, DNA Engine Dyad, Dyad Disciple			•	•	•	•
iCycler	•	•				
iQ, iQ-4, iQ-5, myiQ PCR System	•	•				
MJ Research (Bio-Rad) PTC200, 220, 221, 225, 240	•	•			•	•
MyCycler	•	•				
PTC 200			•	•	•	•
T100	•	•				
Corbett						
Corbett Research (Qiagen) Palm Cycler	•	•				
Corbett Research (Qiagen) Palm Cycler 384					•	•
Eppendorf						
MasterCycler ep 384 / Pro 384					•	•
Mastercycler ep realplex	•	•	•	•		
MasterCycler ep, ep gradient, Pro, Pro S, nexus, nexus gradient, nexus SX1, nexus GSX1	•	•	•	•		
GE						
GE Healthcare Amersham MegaBACE 4000					•	•
GE Healthcare Amersham MegaBACE 500, MegaBACE 1000 mark 2			•	•		
Techne						
Flexigene	•	•	•	•	•	•
Techne Genius, Genius Quad	•	•	•	•	•	•
Techne PCRmax Alpha cyclor 1, 2, 4					•	•
Techne Prime, PrimeG, Prime Elite, Prime Elite Satellite	•	•			•	•
Techne Quantica			•	•	•	•
Techne TC412, TC512			•	•	•	•
Techne TC-PLUS	•	•	•	•		
Touchgene, Touchgene Gradient	•	•	•	•		
Additional Brands						
Apollo ATC-201	•	•				
BIONEER AllInOneCycler	•	•				
G-Storm GS1	•	•	•	•		
HiMedia® Prima-96	•	•				
Idaho Technology LightScanner			•	•	•	•
peqLab peqSTAR 384X						
peqLab peqSTAR 96X	•	•				
Sensoquest Labcycler	•	•			•	•
Takara Dice touch, Gradient	•	•				
VWR Collection UNO96	•	•	•	•		
VWR peqlab peqSTAR 96X			•	•		