


## Serological Pipettes

- Sterile, individually wrapped
- Ascending & descending graduations
- Cotton plugged
- Standardized color coding
- Paper/plastic packaging



MTC Bio serological pipettes are made from virgin polystyrene in a class 100,000 cleanroom facility, under ISO 13485. Available in six volumes (1, 2, 5, 10, 25, 50mL), they are sterilized via E-beam and individually packaged in paper/plastic pouches. All sizes have vivid black, positive and negative graduations and are color coded for easy volume identification.

Cat. No	Description/Volume	Graduations	Color	Qty
P8270-01	Serological Pipette/1mL	1/100	Yellow	500
P8270-02	Serological Pipette/2mL	1/50	Green	500
P8270-05	Serological Pipette/5mL	1/10	Blue	200
P8270-10	Serological Pipette/10mL	1/10	Red	200
P8270-25	Serological Pipette/25mL	2/10	Purple	150
P8270-50	Serological Pipette/50mL	5/10	Black	150

