

UroMAX Series

Urology C-arm Table
by Surgical Tables, Inc



Shown with optional Accessory
package (Stirrups, clamps and drain
loop bag combo)



Standard Accessories:

Hand Control



Foot Control



Patient Straps



Standard Side rails



Easy c-arm access

500 lb weight capacity

2" patient comfort pad

Light weight and easy to move

24" wide top

Made in the USA

Optional Accessories



Drain Loop



Fluid Catchers



Premium Stirrups



Paper Roll



Full Body Extension



Foot Board



Procedural Arm Boards



Surgical Arm boards
(*side rails required)

Most OR accessories can be attached to our side rails

Siderail Options

Split – Standard Style




Straight – Optional Style



Configurations/Specifications


UroMAX-1

Height 28"-38"




UroMAX 3

Height 28"-38"




UroMAX 3

Height 28"-38"




UroMAX-4

Height 30"-40"



UroMAX-5

Height 33.5"-43.5"



500 lb weight limit
Tabletop: 55"L x 24"W
120V/60Hz or 230V/50Hz
5" Locking Swivel Castors
2" Thick Mattress
Lateral Tilt: 17°
Trendelenburg: 12°
Longitudinal Travel 10"
Lateral Travel: 10"
Table Dimensions: 55"L x 26"W