



RC Imaging
your radiology resource

Premium Tabletop Lateral Holder

RC Imaging tabletop holder allows the Radiologic Technologist to position DR panels in a vertical position to achieve those difficult cross-table lateral views.

1.25" channel depth will accommodate most DR panels with or without a panel protector.

The easy glide clamp system makes adjusting vertical height a breeze. This product features non-skid rubber feet to prevent misalignment and ensure quality images.

*Weight 8 Lbs.
Base 11.5"x11.5"
Height 25"
Channel Width 1.25"
Vertical Travel .5"-25.5"*

Part # A.TTP

