



RC Imaging

your radiology resource

Leg Positioning Stand



RC Imaging Leg Positioning Stand converts to 10 telescopic locking positions for effortless vertical adjustments. Lined with pressure sensitive closed cell foam for patient comfort, non-skid rubber feet, and wide steel base designed for added stability. Ideal for obtaining Dunn view images and upright PCL stress views.

- *Cam-Stay technology with 10 locking positions.*
- *Closed cell PSA foam for easy cleaning and sanitation.*
- *Non-skid rubber feet for table grip.*
- *Wide powder coated steel base for patient stability.*
- *Weights 5.6 lbs.*
- *Base Plate: 10 x 10 x .25 in.*
- *Min Height: 13.47 in.*
- *Max Height: 19.00 in.*
- *Leg Holder: 7.5 x 6 x .25 in.*