



July 29, 2021

To Whom It May Concern,

This letter is being sent to confirm the following for Aquasonic[®] 100, Aquasonic[®] Clear, and Scan[®] Ultrasound Transmission Gels:

Aquasonic[®] 100, Aquasonic[®] Clear, and Scan[®] Expiration Date:

Aquasonic[®] 100, Aquasonic[®] Clear, and Scan[®] sold in bottles and tubes have a 5 year shelf life; they expire five years after the manufacture date.

For the original, pre-filled bottles:

Once the bottle has been opened, the gel is expected to maintain its original shelf life of five years, as long as the bottle is stored properly and is not contaminated during use.

For reusable bottles:

Reusable bottles are those that are cleaned and refilled with gel from a larger supply, such as Parker's 5 Liter Aquasonic[®] Ultrasound Gel containers and Scan[®] gallon containers. For these bottles, Parker recommends that the bottles are marked with an expiration date when they are refilled and discarded after one month.

Please let me know if I can be of any further assistance

Best regards,

A handwritten signature in blue ink, appearing to read "Adrina Lalabekova".

Adrina Lalabekova
Quality Assurance Specialist
Parker Laboratories, Inc.
alalabekova@parkerlabs.com