

## MRI/ CT/ X-Ray Testing Summary

Oasis




- MRI Safe
- CT/ Xray Radiopaque

Oasis +



- MRI Safe
- CT/ Xray Radiopaque  
Strongest CT artifact

Azure




- MRI Safe
- Xray Radiopaque
- CT Radiolucent (thin gel top lining visible)

Oasis Elite



- MRI Safe
- CT/ Xray Radiopaque (PU expand foam),
- CT/ Xray Radiolucent (foams invisible- thin gel top lining visible)

ARDS/ Elite Foam



- MRI Safe
- CT/ Xray Radiolucent (Foam invisible- similar to air)

**Radiolucent** substances are materials that allow X-rays or other forms of radiation to pass through relatively easily.

**Radiopaque** substances are materials that do not easily allow X-rays or other forms of radiation to pass through and require higher radiation levels.