

THE 2101KL

DIGITAL PLATFORM SCALE

The all-digital 2101KL is a comprehensive bariatric scale that supports up to 1000 pounds and provides the optimum in quality, reliability, and ease of use. The extra wide “live” handrails allow the patient to hold on for stability while the scale accurately calculates weight.

Designed to meet the needs of doctor’s offices, hospitals, clinics, and acute and long-term care facilities, the scale’s keypad features antimicrobial protection, a display head that swivels 180 °, a large non-skid low-profile platform, and wheels that enable portability.

The scale’s high-contrast color screen provides a user-friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The scale provides easy integration and seamless communication with today’s EMR applications, the 2101KL includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.



THERE’S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it’s not surprising that the 2101KL is your choice for a bariatric platform scale.

120+ LOCAL
FIELD SPECIALISTS

100%
CONNECTIVITY

#1 TRUSTED
BRAND

100 YEARS
OF SERVICE

Health o meter®
Professional

Health o meter® Professional 2101KL

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

1000 lb / 454 kg

Resolution

0.2 lb / 0.1 kg

EMR Connectivity

Via USB, optional Pelstar® wireless technology

Display

1 ½" / 38 mm High-Contrast Color TFT-LCD Screen

Power Source

100-240V adapter (ADPT30) included, 6 D batteries not included

Functions

LB/KG conversion
 LB/KG lock out
 BMI
 Zero
 Tare
 Pre-tare
 Hold/release
 Reweigh
 Recall
 Auto zero
 Auto off
 Time/date
 Variable auto off time
 Audible/mute sound option
 Parent/child weighing function
 Everlock® - permanent weight lock

Platform Size (w x d x h)

26" x 22" x 2 ¾"
 (660 mm x 559 mm x 60 mm)

Product Footprint

31" x 27" x 50"
 (787 mm x 686 mm x 1270 mm)

Product Weight

66 lb / 30 kg

Portability

2 wheels for ease of movement

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
2101KL	Digital platform scale, LB & KG
2101KL-BT	2101KL with Pelstar® wireless technology
2101KG	Digital platform scale, kilograms only
2101KG-BT	2101KG with built-in Pelstar® wireless technology
2101KLHR	2101KL with mechanical height rod included
2101KGHR	2101KG with mechanical height rod included
STROD	Optional mechanical height rod for 2101 models
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
245EHR-2101	Optional digital height rod for 2101 models**

* Not available in all countries

** Compatible with 2101 series scales including a date code with the last two digits of 21 or higher. An upgrade option is available for models including a date code with the last two digits of 17 or higher. For assistance, contact Customer Service at 800-815-6615.

Digital and Mechanical Height Rod Options Available



Health o meter®
 Professional



*Not available in all countries.