



The 2400KL is the perfect scale for locations where occasional wheelchair-bound patients need to be weighed. The scale folds easily and requires minimal room for storage. The scale includes a handle and wheels for easy portability.

Advantages

- EMR Connectivity**
 Scale interfaces with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- Value-Added Features**
 Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option) included at no additional cost
- Calculates Body Mass Index (BMI)**
 Scale functions include the ability to calculate the patient's BMI
- Antimicrobial Keypad**
 The keypad features antimicrobial protection
- Remote Display and Adjustable Tracks**
 Display can be mounted on wall or placed on table. Wheelchair rails expand and adjust from 25 ¾" to 50 ⅛" to accommodate all standard and bariatric wheelchairs
- Easy Mobility**
 Incorporated wheels and casters enable smooth and easy portability, reducing strain on staff
- Convenient Storage**
 Track base folds, can easily be stored in a closet or other small space
- Everlock® (patent pending)**
 Option that allows users to permanently lock the scale in either LB or KG to improve safety



Specifications

- Capacity: 800 lb / 363 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size:
Expands and adjusts from 25 3/4" (654 mm) to 50 1/8" (1273 mm)
- Rail Size (w x d):
6" x 40" / 152 mm x 1017 mm
- Product Footprint, Operating Position with Ramps Fully Extended: (w x h x d):
51 1/4" x 40 1/2" x 3 1/2" / 1302 mm x 1029 mm x 89 mm
- Product Footprint, Folded Position: (w x h x d):
6 3/4" x 40 1/2" x 30" / 171 mm x 1029 mm x 762 mm
- Product Weight: 40 lb / 18 kg
- Wheels: (2) Wheels, (2) Casters
- Display: 1 1/2" / 38 mm High-Contrast Color TFT-LCD Screen
- Connectivity: USB, optional Pelstar[®] wireless technology
- Power Source: 100-240V adapter (ADPT30) included, 6 D batteries (not included)
- Functions: LB/KG Conversion, LB/KG Lock Out, BMI, Zero,Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option, Everlock[®]
- 2 year limited warranty

Model	Description
2400KL	Digital portable wheelchair scale, pounds & kilograms
2400KL-BT	2400KL with Pelstar [®] wireless technology
2400KG	Digital portable wheelchair scale, kilograms only
2400KG-BT	2400KG with Pelstar [®] wireless technology
SS-2400KL	ScaleSurance 2 year extended warranty (2400 models)*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America

* Not available in all countries