

# Health o meter<sup>®</sup>

P r o f e s s i o n a l



✓ **Easy to Use**

Smooth sliding headpiece and telescoping rod allows the user to quickly and easily read the height measurement.

✓ **Large Measurement Range**

The 2" - 90" / 6 cm - 230 cm measuring range provides height readings for a wide range of patients.

✓ **Durable Construction**

Aluminum height rod withstands frequent use. Rigorous laboratory tests showed that the Health o meter<sup>®</sup> Professional STROD lasts 2 ½ times longer than similar models.

### Specifications

- Measuring Range: 2" - 90" / 6 cm - 230 cm
- Graduation: 1/16" / 1mm
- 2 Year Warranty