

THE 2210CART-AC

MULTI-FUNCTION ACUTE CARE CART FOR 2210KL-AM SERIES OF SCALES

- 2 easy grip handles for maneuverability
- Customized plates for secure scale placement
- 2 storage shelves including a built-in 21 ¼"(w) x 13"(d) shelf and a 21 ¼"(w) x 1 ½"(d) hanging shelf
- Four large anti-static casters with locks for easy mobility
- Used to transport neonatal scale or other medical devices and supplies
- Mobile Acute Care cart designed explicitly for the Neonatal/ Pediatric scales 2210KL-AM Series of Scales



2210CART-AC

2210CART-AC Compatible Scales

2210KL-AM
2210KG-AM
2210KL-AM-BT
2210KG-AM-BT
2210KL4-AM
2210KG4-AM
2210KL4-AM-BT
2210KG4-AM-BT
2210KL-AM-C
2210KG-AM-C
2210KL-AM-BT-C
2210KG-AM-BT-C

THE 2210CART

ROLLING SCALE/ MEDICAL DEVICE CART

- The top shelf has a handle for easy grip
- 4 casters (2 casters with locks)
- Spacious storage shelf
- Lockable cabinet
- Used to transport pediatric scales or other medical devices and supplies



2210CART

2210CART Compatible Scales

553KL
553KG
549KL
550KL
2210KL-AM
2210KG-AM
2210KL-AM-BT
2210KG-AM-BT
2210KL4-AM
2210KG4-AM
2210KL4-AM-BT
2210KG4-AM-BT
2210KL-AM-C
2210KG-AM-C
2210KL-AM-BT-C
2210KG-AM-BT-C

Health o meter[®]
Professional

Health o meter[®] Professional 2210CART-AC & 2210CART

weigh**easier**[®]

TECHNICAL SPECIFICATIONS

2210CART-AC MULTI-FUNCTION ACUTE CARE CART FOR 2210KL-AM SERIES OF SCALES

Capacity
77 lb / 35 kg

Product Footprint (w x d x h)
30" x 22" x 38"
762 mm x 559 mm x 965 mm

Storage Shelf (w x d x h)
21 1/4" x 13" x 3"
540 mm x 330 mm x 76 mm

Hanging Shelf (w x d x h)
21 1/4" x 1 1/2"
540 mm x 38 mm

Wheels
4 large anti-static casters with locks

Warranty
2-year limited manufacturer warranty

2210CART ROLLING SCALE / MEDICAL DEVICE CART

Capacity
200 lb / 90 kg

Product Footprint (w x d x h)
26 1/2" x 18" x 34"
673 mm x 457 mm x 867 mm

Cabinet Dimensions (w x d x h)
23" x 16 3/4" x 17"
584 mm x 427 mm x 432 mm

Wheels
4 casters (2 with locks)

Warranty
1-year limited manufacturer warranty

2210CART-AC
Multi-Function Acute Care
Cart for 2210KL-AM Series
of Scales



2210CART
Rolling Scale/ Medical
Device Cart



ORDERING INFORMATION

MODEL	DESCRIPTION
2210CART-AC	Multi-Function Acute Care Cart for 2210KL-AM Series of Scales
2210CART	Rolling Scale/ Medical Device Cart*



*2210CART is a TAA compliant product.

Health o meter®
Professional