

## THE 1100KL-EHR

# DIGITAL PLATFORM SCALE WITH EXTRA WIDE HANDRAILS AND DIGITAL HEIGHT ROD

The all-digital 1100KL-EHR combines the industry-leading bariatric scale that supports up to 1000 pounds and provides the optimum quality, reliability, and ease of use with a state-of-the-art digital height rod. The market-leading height rod provides the most accurate height measurement in its class, with a graduation/resolution of 1mm/0.1" and a full measuring range of 11 ¾" – 79 ½" / 30cm – 202cm. The extra-wide "live" handrails allow the patient to hold on for stability while the scale accurately calculates weight and height.



### Built-In Reset Feature

Proprietary technology eliminates mechanical sensors, reducing failures, increasing accuracy and reliability

Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale's keypad features antimicrobial protection, a display head that swivels 180°, a large non-skid low profile platform, and wheels that enable portability

The scale's high-contrast color screen provides a user-friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare, and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The scale provides easy integration and seamless communication with today's EMR applications. The 1100KL-EHR includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.



## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 1100KL-EHR is the industry leading bariatric scale.

**120+** LOCAL  
FIELD SPECIALISTS

**100%**  
CONNECTIVITY

**#1** TRUSTED  
BRAND

**100** YEARS  
OF SERVICE

**Health o meter®**  
Professional

Health o meter® Professional 1100KL-EHR

Industry leading bariatric scale as reported by Global Industry Analysts, Inc. Patient Scales market survey and Global Healthcare Exchange LLC, (GHX) sales data.

weigh**easier**®

# TECHNICAL SPECIFICATIONS

## Capacity

1000 lb / 454 kg

## Resolution

0.2 lb / 0.1 kg

## EMR Connectivity

Via USB, optional Pelstar® wireless technology

## Display

1 ½" / 38 mm High-Contrast Color TFT-LCD Screen

## Power Source

100-240V adapter (ADPT30) included, 6 D batteries not included

## Digital Height Range

3ft 11.2" – 6ft 7.5" / 120 cm - 202 cm  
Resolution: 1 mm / 0.1"

## Mechanical Height Range

11 ¾" – 79 ½" / 30cm – 202cm  
Graduation: 1 mm / 1/16 inch

## Functions

LB/KG conversion  
LB/KG lock out  
BMI  
Zero  
Tare  
Pre-tare  
Hold/release  
Reweigh  
Recall  
Auto zero  
Auto off  
Time/date  
Variable auto off time  
Audible/mute sound option  
Parent/child weighing function  
Everlock® - permanent weight lock

## Platform Size (w x d x h)

22" x 15 ¾" x 2 ½"  
(559 mm x 400 mm x 64 mm)

## Product Footprint

(height rod extended)  
31 ¾" x 33 ½" x 82 ½"  
(806 mm x 851 mm x 2096 mm)

## Product Footprint

(height rod in home position)  
31 ¾" x 33 ½" x 50 ¼"  
(806 mm x 851 mm x 1276 mm)

## Product Weight

58 lb / 26 kg

## Portability

2 wheels for ease of movement

## Warranty

2-year limited manufacturer warranty

## Industry-Leading Accuracy

The most accurate digital height rod in its class with a graduation/resolution of 1 mm / 0.1"



# ORDERING INFORMATION

MODEL	DESCRIPTION
1100KL-EHR	Digital platform scale with extra wide handrails and digital height rod, LB & KG
1100KL-BT-EHR	1100KL-EHR with Pelstar® wireless technology
SS-1100KL-EHR	ScaleSurance® 2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer

## Height Rod Level to Platform

Bottom of height rod at the platform level, ensuring caregiver confidence and accurate measurement every time

**Health o meter®**  
Professional



\*Not available in all countries.