



8320KL Open

### Advantages

- Motion Detection**  
Scale has the ability to weigh accurately despite movement of the subject on the scale
- Tare Function**  
This convenient function allows the user to automatically subtract the weight of blanket or diaper from infant's weight measurement
- Value-Added Features**  
Operating functions (Tare, Zero Out, Auto Zero, Auto Off) included at no additional cost
- Easy Portability**  
This scale is ideal for traveling nurses or small healthcare facilities with its sturdy construction, built-in carrying handles, easy to clean surface, optional carrying case and ability to fold in half



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
8320KL	4	22.125" x 17.75" x 14.375"/3.27	24.25 LB	10892076002814

**Specifications**

- Capacity: 33 lb / 15 kg
- Resolution: 0.5 oz / 10 g
- 3/4" LCD display
- Functions: LB / KG Conversion, Zero Out / Tare, Auto Zero, Auto Off
- Tray Size: 23-1/2" (w) x 11" (d) x 7/8" (h) / 597mm (w) x 279mm (d) x 22mm (h)
- 2 AAA batteries included
- Optional carrying case (order # 64771)
- 2 Year Limited Warranty



Carrying Case for 8320KL  
Item # 64771