

**PRODUCT INFORMATION**  
P SHIELD GLASSES | LS-SI-BO



The P Shield Glasses is a goggle that provides patient eye protection for Laser and IPL treatments. These P Shield Blockout Glasses are made from soft and pliable TPU, which provides a comfortable fit. In addition, these Patient Blockout Glasses also come with an adjustable strap and a storage case. The Display Wavelength (nm) / OD is 190-10600nm OD 7+. In addition, the P Shield has the following CE Rating:

190-315 D LB8 + IR LB4

>315 1400 D LB6 + IR LB7 + M LB7Y

>1400-11000 DI LB3

**SPECIFICATION**

---

**MATERIAL:** Plastic

**COLOR:** Black

**FEATURES:** Adjustable strap  
and storage case