

PRODUCT INFORMATION

ADVANCE PATIENT BLOCKOUT GLASSES | LS-SS-ADJ-BO



The Patient Blockout Glasses is a stainless steel goggle that provides patient eye protection for Laser and IPL treatments. The Advance Patient Blockout Glasses feature autoclavable glasses and one set of metal rings. In addition, these Patient Blockout Glasses also come with an adjustable strap and a storage case.

OD 7+ 200nm-10600nm
VLT – 0%

SPECIFICATION

COLOR: Black

MATERIAL: Stainless Steel

FEATURES: Adjustable Strap and Storage Case