

# LASER OCULAR SHIELD WITH SUCTION CUP

PRODUCT INFORMATION

COPYRIGHT ©2025. ALL RIGHTS RESERVED



**PRODUCT INFORMATION**  
LASER OCULAR SHIELD WITH SUCTION CUP



The Laser Ocular Shield with Suction Cup is a specialized ocular shield designed to provide optimal eye protection for patients undergoing laser procedures. Crafted with precision from stainless steel, this shield features a non-reflective anterior surface, ensuring a secure and gentle fit. Plus, the addition of a suction cup enables a secure and safe removal process, underscoring its commitment to patient well-being during laser procedures. Available in four sizes (XS, S, M, L) and two orientations (Left, Right), this ocular protector is sold in pairs. With its unique design and attention to detail, the Laser Eye Protector Suction Cup stands as a singular solution and complies with CE EN 93/42/EEC safety standards.

**PRODUCT SPECIFICATION**

---

**COLOR:** Gray

**FEATURES:** Suction Cup

**MATERIAL:** Stainless Steel

**CE Certified**

**AVAILABLE SIZES:**

**Extra Small (22mm x 18mm)**

**Small (26mm x 19mm)**

**Medium (28mm x 21mm)**

**Large (31mm x 22mm)**