



## Advantages

- ✓ **Choice of Units of Measurement**  
 Digital display provides height measurement in centimeters, feet, or inches.
- ✓ **Portability**  
 Lightweight and battery powered, the height rod can be moved anywhere it is needed.
- ✓ **Efficient Design**  
 Measuring infant length is fast and easy with clear digital measurements, large head and foot paddles, and convenient Hold function.

## Specifications

- Measurement Range: 13 ¾" – 31 ½" / 35 -80 cm
- Resolution: 1/16" / 0.1 cm
- Physical Dimensions:
  - Depth: 11 ½" / 290 mm
  - Width (when closed): 24 ¾" / 630 mm
  - Height: 2 ¾" / 70 mm
- Power Source: 2 AAA batteries (included)
- Display: LCD, 5/8" / 16 mm digit height
- Functions: Auto Off, Hold
- 2 year limited warranty

Model	Description
eBABYHR	Digital baby height rod
522EHR	eBABYHR that attaches to 522KL & 522KG
553EHR	eBABYHR that attaches to 553KL & 553KG
550EHRBRACKET	550KL mounting bracket for eBABYHR (must purchase with eBABYHR to mount on 550KL)
SS-eBABYHR	ScaleSurance 2 year extended warranty for eBABYHR*

\* Not available in all countries