



**fitted with**

**work equipment**

**connected load**  
**power consumption**  
**system of protection**  
**class of protection**  
**technology**  
**usage**  
**luminaire head**  
**lamp cover**  
**tubular sections**  
**balance of articulation**  
**joints weight (net)**  
**mains lead**

**fastening**

**lens; dimensions**  
**lens; diopters**  
**add-on lens;**  
**special features**

**3 x 3.5W LED modules, 5000K, 1500 lumens**

**LED convertor**

**120V 50/60 Hz**  
**approx. 13.2 W**  
**IP 20**  
**I**  
**multi-stage switchable**  
**switch**  
**plastic, light grey**  
**screen, milky**  
**steel tube, painted, light grey**  
**spring**  
**approx. 4.5 kg (9.9 lbs)**  
**approx. 3.3 m (11 ft), plug in transformer**

**table clamp (included)**

**162 x 105 (6.38" x 4.13")**  
**3**  
**4 or 4+12 diopters (optional add on) light**  
**output controlled with 2-way-switch**



NEMA 1-15P

**illuminance**

Comment: d=15cm

Em: 4252 lx  
 Emax: 6096 lx  
 Emin: 1967 lx

