

**EASY SLIDE CASSETTE HOLDER
CAS-MB-38479**

PRODUCT INFORMATION

COPYRIGHT ©2025. ALL RIGHTS RESERVED



**PHILLIPS SAFETY
PRODUCTS INC.**

PRODUCT INFORMATION

EASY SLIDE CASSETTE HOLDER | CAS-MB-38479



The Easy Slide Cassette Holder features an easy height adjustment with our unique counter weight design. To adjust the arm height all you simply have to do is slide the Cassette Assembly to your desired distance with minimal effort. This vast improvement over the typical indexed height adjustments style, results in a safe, ergonomic movement for the technician. The Easy Slide also allows for 360 degree rotation of the Cassette with our dual pivoting joints. The Easy Slide Cassette Holder has a heavy mobile base for stability and includes 4 Lockable Casters. This holder ensures effortless mobility and stability in any setting. It accommodates cassettes ranging from 8"x10" to 17"x17", including DR Panel Protectors and Grid Encasements.

Height adjustment from the floor is 6" to 46.5" (Bottom edge of your cassette)

SPECIFICATION

COLOR: Gunmetal

HEIGHT ADJUSTMENT: from the floor 6" to 46.5"
(bottom edge of your cassette)

FEATURES: 4 Lockable Casters, Adjustable Height,
Easy Slide

• Rotates 360 degrees, and tilts +/- 90 degrees allowing cassettes

• Compatible with DR panels and cassettes from 8"x10" to 17"x17" including DR Panel Protectors and Grid Encasements