



### Advantages

**Adjustable Arm & Foot Rest**  
The flip-arm and foot rest allows for easier patient access and patient comfort

**2 Wheels**  
The base of the scale has two wheels that provide easy mobility and allows for convenient placement



| Model | Shipping Dims. (L x W x H) / Cube | Weight (LB) | UPC Code     |
|-------|-----------------------------------|-------------|--------------|
| 445KL | 30.25" x 25.5" x 45" / 20.09      | 81.5 lb     | 892076002046 |



**Specifications**

- Capacity: 440 lb / 200 kg
- Graduation: 1/4 lb / 100 g
- Seat Size: 18 1/2" (w) x 17 1/4" (d) x 15 3/4" (h) / 470 mm (w) x 438 mm (d) x 400 mm (h)
- Wheels: 2
- 2 Year Limited Warranty