



Advantages

- EMR Connectivity**
 Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- Motion Detection**
 Scale has the ability to weigh accurately despite movement of the subject on the scale
- Moveable Arm & Foot Rest**
 The flip-arm and foot rest allows for easier patient access and patient comfort
- Value-Added Features**
 Operating functions (LB/KG Conversion, LB/KG Lock Out, Zero Out/Tare, Auto Off, Reweigh) included at no additional cost



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
2595KL	32.625" x 24.875" x 41.875" / 19.67	100.25 lb	892076002626



Specifications

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Seat Size: 18 1/4"(w) x 15 3/4" (d) x / 464 mm (w) x 400 mm (d)
- Display: 1 1/8" LCD
- Wheels: 2
- Connectivity: USB
- Power Source: 110-240V Adapter (included), 6 C-cell batteries (not included)
- Functions: LB / KG Conversion, LB / KG Lock Out, Zero Out /Tare, Auto Off



2595KL Display