






fitted with	LED - Energy efficiency category A+
Illuminancefield @ 100 cm	Ø 18 cm, (7 in)
illuminance @ 100 cm	100 000 lx
color rendering index (CRI)	> 95
colour temperatur	4300 K
work equipment	electronic converter
mains lead	3 m, (9.84 ft)
luminaire body	polyamide (PA)
weight (net)	20 kg, (44 lb)
Power consumption	40 W
usage	keypad at luminaire head
balance of articulation joints	spring
class of protection	I
Accessory	sterilizable handle
norms	EN 60601-1, EN 60601-2-41
special features	4 dimming levels, Diamond optics, High-quality arm system, sterilisable grip, maintenance free
Fastening	roller stand
System of protection	IP 20 (43)
order number	D15915000 - pure white




light distribution





fitted with	LED - Energy efficiency category A+
Illuminancefield @ 100 cm	Ø 18 cm, (7 in)
illuminance @ 100 cm	100 000 lx
color rendering index (CRI)	>95
colour temperatur	4300 K
work equipment	electronic converter
luminaire body	polyamide (PA)
weight (net)	15.9 kg, (35 lb)
Power consumption	40 W
usage	keypad at luminaire head
balance of articulation joints	spring
class of protection	I
Accessory	sterilizable handle, adapter for suspended ceiling
norms	EN 60601-1, EN 60601-2-41
special features	4 dimming levels, Diamond optics, High-quality arm system, sterilisable grip, maintenance free
Fastening	ceiling mount, when mounting ceiling tube could be cut to length 25-120 cm
System of protection	IP 20 (43)
order number	D16072000 - pure white




light distribution





fitted with	LED - Energy efficiency category A+
Illuminancefield @ 100 cm	Ø 18 cm, (7 in)
illuminance @ 100 cm	100 000 lx
color rendering index (CRI)	>95
colour temperatur	4300 K
work equipment	electronic converter
luminaire body	polyamide (PA)
weight (net)	15.9 kg, (35 lb)
Power consumption	40 W
usage	keypad at luminaire head
balance of articulation joints	spring
class of protection	I
Accessory	sterilizable handle
norms	EN 60601-1, EN 60601-2-41
special features	4 dimming levels, Diamond optics, High-quality arm system, sterilisable grip, maintenance free
Fastening	wall
System of protection	IP 20 (43)
order number	D16077000 - pure white

light distribution





light source
work equipment
connected load
power consumption
color temperature
color rendering index (CRI)
luminous power
illumination area
luminaire body
weight (net)
usage
balance of articulation joints
protection rating

class of protection
standards met
mains lead
accessories
mounting
special features

medical device class
battery life cycle up to 70%
battery type

LED module
electronic plug-in transformer
100-240 V, 50/60 Hz
40 W
4300 K
>96
100 000 lux / 1.0 m
18cm (7.1 in) @ 1.0m
polyamide (PA)
70.5 lbs (32 kgs)
keypad at luminaire head
spring (weight compensation)
IP 43 (luminaire head)
IP 20 (UPSmed pack)
I
EN 60601-1, EN 60601-2-41
3 m cord (9.8 ft)
sterilizable hand, handle cover
roller stand
diamond optics, sterilizable handle,
4 dimming levels
I according to MDR
1,200 cycles
NP17-12I industrial VRLA battery

light distribution

