



TRIANGO 60
THE COMPACT ONE THAT HAS WHAT IT TAKES



TRIANGO 60 – your new treatment & minor surgery light



TRIANGO 60 C

ceiling version



TRIANGO 60 W

wall mount



TRIANGO 60 F

roller stand



TRIANGO 60 FB*

roller stand



| Mounting Option | Arm Length (in) |
|-----------------|-----------------|
| ceiling mount | 68 |
| wall mount | 78 |
| rolling stand | 36 |

| TRIANGO 60 | technical data |
|---------------------------------------|--|
| Fitted with | 1 × LED-Modul |
| Luminous power ** | 60000 Lux / 1.0 m |
| Color temperature ** | fixed value 4700 K |
| TRIANGO 60-1 F | changeable value |
| TRIANGO 60-3 F | 3700 / 4300 / 4700 K |
| Lighting field 1.0 m** | ∅ 18 cm (7 in) ∅ 18/23/25 cm (7/9/10 in) focusable option |
| Depth of illumination ** | 140 cm (55 in) |
| Color rendering index | CRI > 95 |
| Operation equipment | Electronic converter |
| Voltage | 100-240 V; 50/60 Hz |
| Luminaire body | Polyamide (PA) |
| Power consumption | < 30 W |
| Usage | Touch panel on cardanic system |
| Technology | Single-pole secondary |
| Balance of articulation joints | Spring (weight compensation) |
| Class of protection | I |
| Protection rating | IP 43 (illumination head) IP 20 (UPSmed Pack) |

| Special features | |
|--------------------------------|---|
| TRIANGO 60-1 | <ul style="list-style-type: none"> • Continuously dimmable |
| TRIANGO 60-3 | <ul style="list-style-type: none"> • Continuously dimmable • 3 color temperatures |
| TRIANGO Fokus 60-1 | <ul style="list-style-type: none"> • Continuously dimmable • Focusable |
| TRIANGO 60-1 F Battery* | <ul style="list-style-type: none"> • Continuously dimmable • UPSmed Pack |
| Mounting | <ul style="list-style-type: none"> • Ceiling mount • Wall mount • Roller stand |
| Cord length | |
| TRIANGO 60-1 C | Direct power supply |
| TRIANGO 60-1 W | Direct power supply |
| TRIANGO 60-1 F | 3 m (9.8 ft) |
| TRIANGO 60-1 F Battery* | 3 m (9.8 ft) |
| Weight (net) | |
| TRIANGO 60-1 C | 13 kg (28.7 lb) |
| TRIANGO 60-1 W | 12.5 kg (27.6 lb) |
| TRIANGO 60-1 F | 19 kg (41.9 lb) |
| TRIANGO 60-1 F Battery* | 31 kg (68.3 lb) |
| Accessories | <ul style="list-style-type: none"> • Sterilizable handle • Disposable cover for handle |

* Coming soon

** -10% / +20% tolerance

Standards: ANSI/AAMI ES60601-1:2005 / EN/IEC 60601-1 / EN/IEC 60601-2-41 / EN/IEC 60601-1-2 / CAN/CSA-22.2 NO. 60601-1 / CAN/CSA-C22.2 NO 60601-2-41 / CAN/CSA-C22.2 NO. 60601-1-6

Compliance: RoHS

Conformity: CE

Listed: MET/cMET

Versatile

versions with focusable or color changeable illumination field

Homogeneous light field

thanks to unique diamond optic

Easy cleaning

streamlined design and the removable handle

Hygiene made simple

good design does without unnecessary details and attaches importance to closed building structures

Can be used everywhere

TRIANGO 60 is suitable for the versatile daily routine in family practices and clinics as well as in gynecology, dermatology or veterinary medicine

**Ergonomic handling**

easy alignment for right- and left-handed

**Intuitive operation**

thanks to an intuitive touch panel on the cardanic, the illumination from 60,000 lux at 1 m can be infinitely adjusted

Order number: TRIANGO 60-1

| | |
|--------------------------------|-----------|
| TRIANGO 60-1 C | D16253000 |
| TRIANGO Fokus 60-1 C | D16255000 |
| TRIANGO 60-1 W | D16274000 |
| TRIANGO Fokus 60-1 W | D16276000 |
| TRIANGO 60-1 F | D16265000 |
| TRIANGO Fokus 60-1 F | D16267000 |
| TRIANGO 60-1 F Battery* | D16280000 |

Order number: TRIANGO 60-3

| | |
|-----------------------|-----------|
| TRIANGO 60-3 C | D16254000 |
| TRIANGO 60-3 W | D16275000 |
| TRIANGO 60-3 F | D16266000 |

5 YEARS WARRANTY
for all electrical components


MDR CERTIFIED
European Medical Device Regulation


SWISS Engineering


MADE IN SWISS

