




## Advantages

- ✓ **Adjustable Unit Measurement Bar**  
The mechanical rotating poise bars allow for easy reading of weight in lb or kg
- ✓ **Engraved Numerals**  
Engraved numerals on the poise bars can withstand frequent use, providing long-lasting readability
- ✓ **Convenient Weighing Seat**  
Plastic, easy-to-clean seat limits infant movement while accurately calculating weight measurement

**Health o meter**<sup>®</sup>  
P r o f e s s i o n a l



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
1524KL	25" x 20.625" x 18.625" / 5.56	39.5 lb	892076002183

**Specifications**

- Capacity: 130 lb / 65 kg
- Graduation: 1 oz / 20 g
- Seat Size: 24 1/2"(w) x 14 1/4" (d) x 18 1/2" (h) / 622 mm (w) x 362 mm (d) x 470 mm (h)
- Easy to read rotation poise bars for dual read in lb and kg
- Engraved numerals
- 2 Year Limited Warranty

**Health o meter**<sup>®</sup>  
Professional  
[www.homscales.com](http://www.homscales.com)