



This scale is ideal for infants and with its mechanical height rod practitioners can measure both weight and height simultaneously.

### Advantages

- ✓ **EMR Connectivity**  
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ✓ **Motion Detection**  
Scale has the ability to weigh accurately despite movement of the subject on the scale
- ✓ **Scale Value-Added Features**  
Operating functions (LB/KG Conversion, KG Lock Out, Zero Out/Tare, Hold/Release, Auto Off) included at no additional cost
- ✓ **Mechanical Height Rod Included**  
Easy to read height rod allows user to measure both weight and height simultaneously



### Scale Specifications

- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- Tray Size (w x d x h):  
24 7/8" x 13" x 1 3/8"  
632 mm x 330 mm x 35 mm
- Display: 1" / 25 mm LCD
- Power Source: 6 AA batteries included,  
optional power adapter  
(order # ADPT30)
- Connectivity: USB, optional Pelstar<sup>®</sup>  
wireless technology
- Functions: Zero, Tare, Hold / Release,  
Auto Off
- 2 year limited warranty

### Height Rod Specifications

- Range of measurement:  
7 cm – 100 cm / 3" – 39"
- Division of measurement:  
1 mm / 1/16"

| Model        | Description   |
|--------------|---|
| 553KL-HR     | Digital pediatric scale (pounds & kilograms) with mechanical height rod |
| 553KG-HR     | Digital pediatric scale (kilograms only) with mechanical height rod     |
| 2210CART     | Rolling cart for 553KL-HR & 553KG-HR                                    |
| ADPT30       | Optional 100-240V adapter for scale                                     |
| PLUGSET30    | Plug adapters for use outside North America for scale                   |
| PELSTARONE-P | Wireless thermal printer for scale                                      |