

SHIP FLAT BARRIER

PRODUCT INFORMATION



PHILLIPS SAFETY
PRODUCTS INC.

PRODUCT INFORMATION
SHIP FLAT BARRIER



The Ship Flat Barrier is a lightweight lead mobile barrier made for mobility, maximum protection, and global shipping. This barrier features a lightweight, super-high-strength aluminum frame and hospital-grade casters for easy mobility, with 4 brakes included. Weighing 135 lbs, it ships conveniently via UPS.

The barrier also features a lead acrylic window measuring 30 inches wide by 24 inches high with .50mm L.E. on the top and bottom, and a lower opaque panel measuring 30 inches wide by 42 inches high with .5mm L.E. Additionally, assembly is easy with just one box, one wrench, and four bolts (all included).

SPECIFICATION

BARRIER MATERIAL: Aluminium

FEATURES: Lead window, lower opaque panel,
Hospital-grade, four with brakes

OVERALL SIZE: 72" H x 34½" W x 25" D

PROTECTION LEVEL: 0.50mm PbEq top and bottom

WINDOW MATERIAL: Acrylic

WINDOW SIZE: 30" x 24"

PANEL SIZE: 30" x 42"

SHIPS UPS: 135 lbs