

Health o meter[®]

P r o f e s s i o n a l



Advantages

- Weighs Infants and Toddlers**
Dual-purpose scale can be used with the infant weighing tray or toddler step-on platform.
- Weight Hold**
Weight value is held on the display so caregiver can attend to patient before recording the weight.
- Imperial and Metric Measurement**
Offers users the opportunity to weigh patients in whichever unit is required.
- Built-In Connection Port for eBABYHR**
Digital baby height rod simply snaps onto the scale for a comprehensive weighing and measuring device.
- Removable Tray**
The removable tray makes it easy to clean the device.
- Tare Function**
Conveniently subtracts the weight of a blanket or diaper from the patient's weight.
- Easy to Use**
Lightweight and portable with simple functionality and easy to read digital display.



550KL



Specifications

- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- Tray Size (w x d x h):
23 1/2" x 8" x 2"
597 mm x 203 mm x 51 mm
- Platform Size: (w x d x h):
12 1/2" x 8 5/8" x 1 1/4"
318 mm x 219 mm x 32 mm
- Product Footprint Tray (w x d x h):
23 1/2" x 13 1/4" x 3 5/8"
597 mm x 337 mm x 92 mm
- Product Footprint Platform: (w x d x h):
13 1/2" x 11 1/2" x 1 3/4"
343 mm x 292 mm x 44 mm
- Display: 1/2" / 13 mm digit height LCD
- Power Source: 9V battery included
- Functions: Auto Off, Hold, Tare, Unit
- 2 Year Limited Warranty

| Model | Description |
|---------------|---|
| SS-550KL | ScaleSurance 2 year extended warranty* |
| eBABYHR | Digital baby height rod for 550KL (must purchase with 550EHRBRACKET to mount on scale) |
| 550EHRBRACKET | 550KL mounting bracket for eBABYHR (must purchase with eBABYHR to mount on scale) |
| 553CASE | Optional carrying case for 550KL |
| 2210CART | Optional rolling cart for 550KL |

* Not available in all countries