


NOTES:

MATERIAL:
STEEL WITH 2.0mm LEAD
EQUIVALENCE


PAD PRINT COLOR:
2.0mm LEAD
EQUIVALENCE GLASS

FRAME COLOR:
TAN


NOT TO SCALE



A
1 BARRIER ELEVATION



B
1 GLASS ELEVATION



C
1 SECTION AT FRAME

DWG No.	LB-3672	Rev.
P/N:	LB-3672	
DESC.	Mobile Leaded Barrier 36 in. x 72 in Window	