


10" DEEP OPENING AT TOP



14 3/4" INSIDE DIMENSION
(LEFT TO RIGHT)

NOTES

MATERIAL: 22 GA. (.030") 430 STAINLESS STEEL
HOLDS GLOVES, MASKS AND GOWNS
OUTSIDE DIMENSIONS: 26"H x 21"W x 4 1/2"D.
OPTIONAL 307001 - DOOR HANGER BRACKETS
4 KEYHOLE SLOTS PROVIDING FOR MOUNTING. USE UP TO #10 HARDWARE
(NOT INCLUDED). 16" SQUARE PATTERN; 4 7/8" FROM TOP EDGE.

CAT. NO.: 307000		
DESCRIPTION: ISOLATION STATION		
SCALE: NONE	DRAWING NUMBER:	DRAWN BY: R.Z.
CODE: 071601	307000 SPECS.	REVI 6D: