

Improved!

Stainless Steel Support

Stainless steel support jacks from BrandTech Scientific are rugged supports that are excellent for a wide variety of critical laboratory applications. Constructed entirely from stainless steel with POM rollers and PF31 phenolic knob, these support jacks are exceptionally stable and well-suited for corrosive environments such as fume hoods, water baths, marine and industrial labs, and more.

- **Supports a Wide Variety of Items:** Stainless steel support jacks are available in a range of sizes and load capacities up to 176 pounds for most lab applications.
- **Improved Lift Mechanism:** Provides increased stability and torsional rigidity.
- **Helps Reduce Contamination:** Support jacks feature smooth surfaces that are easy to clean and sanitize. All support jacks are fully autoclavable at 121°C for sterile applications.
- **Reaches New Heights:** These stainless steel support jacks are now also available in Maxi-Lift models. Maxi-Lift stainless steel support jacks have a maximum height that is approximately 30% higher than traditional support jacks for special applications.
- **Offers a Wide Variety of Options:** BrandTech Scientific offers support jacks manufactured with anodized aluminum and powder-coated aluminum decks (on stainless steel lift mechanisms) for less critical applications.
- **Two-year warranty on all support jacks**



Rugged Laboratory Supports for Critical Environments.

Jacks

Optional Handle

An optional ratcheting handle is available to simplify use of BrandTech's larger support jacks.

This handle simply slides onto the handwheel of 30 x 30cm and 40 x 40cm support jacks to greatly reduce the effort required for raising and lowering the instrument.



Also Available

Support jacks with epoxy powder-coated aluminum and anodized aluminum decks, with stainless steel lift mechanisms.

Epoxy Powder Coat



Anodized Aluminum



ORDERING INFORMATION

Description Plate, L x W, cm	Minimum Height, cm	Maximum Height, cm	Maximum Static Load, kg	Cat. No.
Stainless Steel Support Jacks				
10 x 10	5.5	12	10	B11115
16 x 13	6.0	27.5	20	B11120
20 x 20	6.0	27.5	20	B11130
24 x 24	6.0	27.5	20	B11140
30 x 30	13.0	47	80	B11180
40 x 40	13.0	47	80	B11190

Maxi-Lift Stainless Steel Support Jacks

16 x 13	7.5	40	20	B11122
20 x 20	7.5	40	20	B11132

Anodized Aluminum Deck Support Jacks

10 x 10	5.5	12	10	B11016
16 x 13	6.0	27.5	20	B11021
20 x 20	6.0	27.5	20	B11031
24 x 24	6.0	27.5	20	B11041
30 x 30	13.0	47	80	B11081
40 x 40	13.0	47	80	B11091

Powder-coated Aluminum Deck Support Jacks

10 x 10	5.5	12	10	B11015
16 x 13	6.0	27.5	20	B11020
20 x 20	6.0	27.5	20	B11030
24 x 24	6.0	27.5	20	B11040
30 x 30	13.0	47	80	B11080
40 x 40	13.0	47	80	B11090

Accessories

Optional Ratchet Wrench

for B11180, B11190, B11081, B11091, B11080, B11090

B11089

1cm = 0.394 inches, 1kg = 2.21 pounds