

Sample Cups for Clinical Analyzers

BRAND sample cups are made from virgin polystyrene and are designed to fit many common analyzers. Caps and lids are made of polyethylene for a snug fit. Cups and caps meet CE-IVD Directive 98/79 (except 722060).



[a] Sample Cups for Technicon® Analyzers



[b] Push-on Caps to fit cup #115015



[c] Sample Cups for Coulter Counter® with Lids



[d] Sample Cups with Snap Caps

ORDERING INFORMATION

	Capacity mL	Upper Diameter mm	Lower Diameter mm	Height mm	Qty./Pack	Cat. No.
Sample Cups for Technicon® Analyzers [a]	1.5	15	12.2	22.7	12000	115015
	2	15.5	12.5	23.5	10000	115016
	4	17	13.3	38	6000	115017
Push-on caps to fit cup #115015 [b]					1000	115020
Sample Cups for Coulter Counter® with Lids [c]	20	32	n/a	56	1000	722055
Sample Cups with Snap Caps [d]	12	22	22	38	1000	722060