

# DETECTO 750 Clinical Weight Indicator

## SERIAL CONNECTIVITY STANDARD

Serial ports on all model 750 indicators used with DETECTO medical scales.



### Serial I/O Connector

RS232 serial port allows for electronic medical recordkeeping



### Scale Connectivity for Patient Records

DETECTO's model 750 digital weight indicators provide RS232 serial connectivity as a standard feature. This allows any DETECTO medical scale utilizing the 750 indicator to have EMR/EHR-ready serial output for interfacing with a computer, monitor, or printer for efficient data transfer in electronic medical recordkeeping.




- Body Mass Index calculation
- Extra-large, 0.75 in/19 mm high LCD digits
- 6 AA batteries (not included) or optional AC adapter
- Lock/Release key for viewing weight after the patient steps off the scale
- Units conversion for lb or kg
- Durable steel construction for years of reliable service



# 750 INDICATOR

STANDARD ON THESE DETECTO SCALE MODELS:



**6127**  
Waist-High Scale



**6129**  
Waist-High Scale  
With Height Rod



**6337**  
Physician Scale



**6339**  
Physician Scale  
With Height Rod



**6475**  
Chair Scale



**6875**  
Chair Scale



**6880**  
Chair Scale



**6800**  
Bariatric Scale



**6855**  
Bariatric Scale



**6856**  
Bariatric Scale



**6876**  
Bariatric  
Flip Seat Scale



**6868**  
Bariatric  
Flip Seat Scale



**IBFL500**  
In-Bed Scale



**IB800**  
In-Bed Scale



**8432CH**  
Pediatric Scale



**8435**  
Pediatric Scale



**DR400-750**  
Platform Scale



**GP400-750**  
Platform Scale



**VET50**  
Veterinary Scale



**VET400**  
Veterinary Scale