



## 1.2 mL capacity – short stem




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
138040	Non-sterile	500 pcs/dispenser box, 10 boxes/case	5000	<ul style="list-style-type: none"> <li>Capacity: 1.2 mL</li> <li>Total length: 65 mm (2.56")</li> <li>Drops/mL: 25</li> <li>Stem diameter: 2.5 mm (0.1")</li> <li>Bulb draw: 0.9 mL</li> <li>Tip diameter: 2.5 mm (0.1")</li> </ul>
	Non-sterile	500 pcs/dispenser box	500	

## 4 mL capacity – short stem – large bulb




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
136020	Non-sterile	500 pcs/dispenser box, 10 boxes/case	5000	<ul style="list-style-type: none"> <li>Capacity: 4 mL</li> <li>Total length: 85 mm (3.4")</li> <li>Drops/mL: 24</li> <li>Stem diameter: 2.5 mm (0.1")</li> <li>Bulb draw: 3 mL</li> <li>Tip diameter: 2.5 mm (0.1")</li> </ul>
	Non-sterile	500 pcs/dispenser box	500	

## 4.5 mL capacity – standard



Catalog No.	Sterility	Packaging Information	Unit	Technical Data
136030	Non-sterile	500 pcs/dispenser box, 10 boxes/cases	5000	<ul style="list-style-type: none"> <li>Capacity: 4.5 mL</li> <li>Total length: 155 mm (6.1")</li> <li>Drops/mL: 25</li> <li>Stem diameter: 2.5 mm (0.1")</li> <li>Bulb draw: 3.2 mL</li> <li>Tip diameter: 2.5 mm (0.1")</li> </ul>
	Non-sterile	500 pcs/dispenser box	500	


## 6.5 mL capacity



Designed for use with spectrophotometer cuvettes

Catalog No.	Sterility	Packaging Information	Unit	Technical Data
139020	Non-sterile	500 pcs/dispenser box, 10 boxes/case	5000	<ul style="list-style-type: none"> <li>Capacity: 6.5 mL</li> <li>Total length: 155 mm (6.1")</li> <li>Drops/mL: 28</li> <li>Stem diameter: 7.8 mm (0.31")</li> <li>Bulb draw: 3.2 mL</li> <li>Tip diameter: 2.1 mm (0.08")</li> </ul>
	Non-sterile	500 pcs/dispenser box	500	

## 15 mL capacity – large bulb



Catalog No.	Sterility	Packaging Information	Unit	Technical Data
139040	Non-sterile	250 pcs/dispenser box, 10 boxes/case	2500	<ul style="list-style-type: none"> <li>Capacity: 15 mL</li> <li>Total length: 155 mm (6.1")</li> <li>Drops/mL: 19</li> <li>Stem diameter: 3 mm (0.12")</li> <li>Bulb draw: 8 mL</li> <li>Tip diameter: 3 mm (0.12")</li> </ul>
	Non-sterile	250 pcs/dispenser box	250	