



4 mL capacity – 5 inch




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
138060	Non-sterile	500 pcs/dispenser box, 10 boxes/case	5000	<ul style="list-style-type: none"> Capacity: 4 mL Total length: 130 mm (5.1") Drops/mL: 25 Stem diameter: 6.6 mm (0.26") Bulb draw: 2 mL Tip diameter: 2.5 mm (0.1")
	Non-sterile	500 pcs/dispenser box	500	

5 mL capacity – 6 inch




Catalog No.	Sterility	Packaging Information	Unit	Technical Data
138050	Non-sterile	500 pcs/dispenser box, 10 boxes/case	5000	<ul style="list-style-type: none"> Capacity: 5 mL Total length: 155 mm (6.1") Drops/mL: 23 Stem diameter: 6.8 mm (0.27") Bulb draw: 2 mL Tip diameter: 2.6 mm (0.1")
	Non-sterile	500 pcs/dispenser box	500	

5 mL capacity – 6 inch – graduated to 2 mL



Catalog No.	Sterility	Packaging Information	Unit	Technical Data
137040	Non-sterile	500 pcs/dispenser box, 10 boxes/case	5000	<ul style="list-style-type: none"> Capacity: 5 mL Total length: 155 mm (6.1") Drops/mL: 23 Stem diameter: 6.5 mm (0.25") Bulb draw: 1.8 mL Tip diameter: 2.8 mm (0.11") Graduated up to: 2 mL Graduation increments: 0.5 mL
	Non-sterile	500 pcs/dispenser box	500	

5 mL capacity – graduated to 2 mL



Catalog No.	Sterility	Packaging Information	Unit	Technical Data
137050	Non-sterile	500 pcs/dispenser box, 10 boxes/case	5000	<ul style="list-style-type: none"> Capacity: 5 mL Total length: 155 mm (6.1") Drops/mL: 23 Stem diameter: 6.8 mm (0.27") Bulb draw: 2 mL Tip diameter: 2.6 mm (0.1") Graduated up to: 2 mL Graduation increments: 0.5 mL
	Non-sterile	500 pcs/dispenser box	500	