

Disposable LUCA® & Receptacle System

INSTRUCTIONS FOR USE



1

Put roll of LUCA® in Receptacle feeding from the top



2

Pull LUCA® until it reaches Perforation



3

Press down gently on Receptacle lid with one hand



4

Tear gently across perforation

 **CAUTION: DO NOT USE LIQUIDS EXCEEDING 37°C**

Please keep a copy of IFU - this must be retained for future reference.

CMCU 10001A - PathAguard® LUCA®

CMCU 10002A - PathAguard® LUCA® Large

CMCD 10002A - PathAguard® LUCA® Receptacle

