

Phillips Safety Products  
123 Lincoln Boulevard  
Middlesex,  
NJ 08846, USA  
01.04.2019

### **Transfer of certification from SGS United Kingdom Ltd (NB 0120) to SGS Fimko Oy (NB 0598)**

Following your request, I confirm that the Article 10/Module B certificates listed in Annex A to this letter have been formally transferred to SGS Fimko Oy, Helsinki, who will act as the European Notified Body. The Notified Body Number for SGS Fimko Oy, listed on the European Commission's NANDO database is 0598.

The Certificate Number has not changed, and retains the prefix associated with SGS United Kingdom Ltd. who originally issued the documents.

Although we cannot mandate any specific text for you to use in your Declaration of Conformity, we suggest that when listing the certificate number in the declaration, you indicate both the old UK and the new EU-27 Notified Body's details / Identification numbers.

Yours sincerely



K R Shepherd  
(CRS) PPE Product Manager, SGS United Kingdom Ltd  
Authorised signatory on behalf of SGS Fimko Oy

**Annex A**

<b>Article 10 certificate numbers</b>	<b>Module B certificate numbers</b>	<b>Product Description/Code</b>
GB17/873435		X-Ray radiation filtering glasses
GB14/92296		X-Ray radiation filtering glasses
GB14/92297		X-Ray radiation filtering glasses
GB1492298		X-Ray radiation filtering glasses
GB17/872574		X-Ray radiation filtering glasses
GB17/872578		X-Ray radiation filtering glasses
GB17/872581		X-Ray radiation filtering glasses
GB14/92293		X-Ray radiation filtering glasses
GB14/92294		X-Ray radiation filtering glasses

# Phillips Safety Products, Inc.

123 Lincoln Boulevard,  
Middlesex,  
New Jersey 08846,  
USA.

It is certified that the manufacturer's technical file and the PPE product detailed on page 2 have been assessed and found to be in accordance with

## EC Directive 89/686/EEC EC Type-examination

This certificate is valid from 06 November 2017 until 06 November 2022  
Issue 1. Certified since 06 November 2017  
Previously certified on certificates GB14/92290, GB14/92291 &  
GB14/92292


Authorised by

SGS United Kingdom Limited, Notified Body 0120  
Unit 202B Worle Parkway, Weston-super-Mare, BS22 6WA, United Kingdom  
t +44 (0)1934 522917f +44 (0)1934 522137  
e-mail [sgsprodcert@sgs.com](mailto:sgsprodcert@sgs.com) [www.sgs.com](http://www.sgs.com)

SGS EC 03 0311

This certificate remains the property of SGS United Kingdom Ltd to whom it must be returned on request.

Page 1 of 2





## Phillips Safety Products, Inc.

### EC Directive 89/686/EEC

#### EC Type-examination

Issue 1.

#### PPE Product

- PS (logo) RG series frames, fitted with X-ray radiation filtering lead glass lenses, with a lead equivalence of 0.75mm Pb at 150 kV;
- PS (logo) RG-75SS series frames, fitted with X-ray radiation filtering lead glass lenses, and side shields, both with a lead equivalence of 0.75mm Pb at 150 kV;
- PS (logo) RG-50SS series frames, fitted with X-ray radiation filtering lead glass lenses, with a lead equivalence of 0.75mm Pb at 150 kV, and side shields with a lead equivalence of 0.50mm Pb at 150 kV;
- PS (logo) RG50 series frames, fitted with X-ray radiation filtering lead glass lenses, with a lead equivalence of 0.50mm Pb at 150 kV;
- PS (logo) RG50-50SS series frames, fitted with X-ray radiation filtering lead glass lenses, and side shields, both with a lead equivalence of 0.50mm Pb at 150 kV.

The frames are available in a variety of designs and colours, which have no bearing on the protection offered.

The following lens options are available for each frame: Plano (no suffix code), Single vision prescription (suffix: -SV), Bifocal (suffix: -FT-28) & Progressive (-PROG).

Lens coatings available are: Fog-Free (-FF) & Anti-reflective (-AR).

It is certified that the manufacturer's technical file and the above mentioned PPE have been assessed and found to be in accordance with the requirements of Council Directive 89/686/EEC. When examined the PPE satisfied the relevant requirements in the manufacturer's technical specification<sup>1</sup> based on (based on EN 61331-3:2014(Protective devices against diagnostic medical X-radiation – Part 3: Protective clothing and protective devices for gonads.) and on testing to EN 61331-1:2014 plus guidance notes) which had been previously assessed and deemed to meet the basic health and safety requirements (Annex II) of the Directive.

Continued...

# Phillips Safety Products, Inc.

## EC Directive 89/686/EEC EC Type-examination

Issue 1.

The manufacturer's information notice (Appendix 5 User Guide Insert) in English, has been inspected and found to have addressed all of the relevant requirements of the Directive (Annex II, 1.4). However, the detailed content of the information notice is the responsibility of the manufacturer, as are translations into other languages.

This certificate is issued on the strict condition that appropriate checks on manufactured PPE, as detailed in Article 11 of the Directive are implemented and maintained while the model is in production.

Certification is based on technical file Radiation Reducing Eyewear dated 4/12/2017  
SGS Reference Number UK/CRS 222428.

<sup>1</sup> The validity of this certificate may have to be reviewed in the light of new or amended harmonized standards.

This certificate remains the property of SGS United Kingdom Ltd to whom it must be returned on request.